



PLEASE READ INSTRUCTIONS CAREFULLY BEFORE USE

- FDA Approved
- 99% Accurate
- Independently Tested
- Scientifically Proven

For more information visit www.testsure.com


ONE STEP STRIP STYLE HCG URINE PREGNANCY TEST

INTRODUCTION

One Step Strip Style HCG Urine Pregnancy Test is a test kit for the determination of HCG (Human Chorionic Gonadotropin) in urine specimens. This test kit is used to obtain a visual qualitative result for early detection of pregnancy. The test sensitivity is 25-30mIU/mL.


TEST PROCEDURE

1. To begin testing, open the sealed pouch by tearing along the notch. Remove the test form the pouch.
2. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse past the MAX (maximum) line. Take the strip out after 5-7 seconds and lay the strip flat on a clean, dry, nonabsorbent surface (e.g., mouth of the urine container).
3. Wait for colored lines to appear. Depending on the concentration of HCG in the test specimen, positive results may be observed in as short as 40 seconds. However, to confirm negative results, the complete reaction time (5 minutes) is required. Do not read results after 10 minutes.



INTERPRETATION OF RESULTS

- Negative:** Only one color line appears on the control region. No apparent line on the test region means you are not pregnant.
- Positive:** Distinct color line appears on the control and test regions. Both test line and control line indicate that you are pregnant. Color intensity of the test lines may vary. Different stages of pregnancy have different concentrations of HCG hormone.
- Invalid:** No visible line at all. Repeat test with a new test kit.



STORAGE AND STABILITY

The test kits can be stored at room temperature (18 to 30C) in the sealed pouch to the date of expiration. The test kits should be kept away from direct sunlight, moisture and heat.

PRECAUTION

1. For in vitro diagnostic use only.
2. Do not use test kit after the expiration date.
3. The test device should not be reused.

ONE STEP STRIP STYLE LH URINE OVULATION TEST

INTRODUCTION

The One Step LH Urine Ovulation Test is fast and easy to use. It is a qualitative test that can predict when there is a LH (luteinizing hormone) surge, and in turn, when you are likely to ovulate.

WHEN TO BEGIN TESTING

First, you must determine the length of your menstrual cycle. This is the number of days from the first day of your menstrual bleeding to the day before your next bleeding begins. Please refer to the chart to determine when you should start testing. If your cycle is shorter than 21 days or longer than 38 days, consults with you doctor. If you do not know your cycle length, you may begin the test 11 days after your first period since the average cycle length is 28 days. Perform 1 test each day over a 5 day period, or until the LH surge has been detected.


Your Cycle Length	Start To Test On
21 days	Day 6
22 days	Day 6
23 days	Day 7
24 days	Day 7
25 days	Day 8
26 days	Day 9
27 days	Day 10
28 days	Day 11
29 days	Day 12
30 days	Day 13

URINE COLLECTION

1. Do not use the early morning urine samples as LH is synthesized in your body during the morning. It will not show up in your urine until later in the day.
2. The best time to collect your urine is between 10am and 8pm.
3. Collect urine at about the same time each day.
4. Reduce liquid intake about 2 hours before collecting your urine as a diluted urine sample can prevent the test from detecting the LH surge.


TEST PROCEDURE

1. Allow the sealed LH test pouch, urine sample, and control solution to reach room temperature (18-30 degrees Celsius)
2. Remove the LH test strip from the sealed pouch.
3. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse past the MAX (marker line). Take the strip out after 5 seconds and lay the strip flat on a clean, dry, non-absorbent surface (such as the mouth of the urine container)
4. Wait for colored lines to appear. Depending on the concentration of LH in the test specimen, positive results may be observed in as little as 40 seconds. However, to confirm negative results, the complete reaction time of 10 minutes is required. Do not read results after more than 30 minutes.



INTERPRETATION OF RESULTS

- Negative:** Only one color line appears on the control region, or the test line appears but is lighter than the control line. This means there is no LH surge.
- Positive:** If two color lines are visible, and the test line is equal to or darker than the control line, one will probably ovulate in the next 24-48 hours. If trying to get pregnant, the best time to have intercourse is after 24 but before 48 hours.
- Invalid:** No visible line at all. Repeat with a new test kit.



STORAGE AND STABILITY

The test kit can be stored at temperature between 2 to 30 degrees Celsius in the sealed pouch to the date of expiration. The test kit should be kept away from direct sunlight, moisture and heat.